

Thwing-Albert

INSTRUMENT COMPANY

ASP-2.2 Semi-Auto Sheet Press

The ASP-2.2 is designed to press more than the suggested number of handsheets that the standard press can handle. A large space has been provided between the upper and lower press platens to allow you to press up to 30 handsheets at a time.

Programmed to meet both CPPA and TAPPI standards, the microprocessor has been programmed to take 5.5 min. for the first pressing and 2.5 min for the second pressing.

The unit has been designed to facilitate easy loading of the handsheets and blotters for pressing. The bottom press platen slides forward and allows the use of the loading template to ensure that the handsheets are centered in the press.

As a safety feature, the act of lowering the top press platen requires the operator to press the two "Down" buttons at the same time for a minimum of 5 seconds and until the top platen touches the top of the handsheet stack. This ensures that the operator will not catch a hand or finger in the press section.

Features

- Complies with CPPA and TAPPI
- 30 handsheets at a time

Specifications

- Physical Dimensions
35 inch H x 16 inch W x 22 inch D
(889 mm x 406.4 mm x 558.8 mm)
- Net Weight: 80 lb (36.4 kg)



▲ ASP-2.2 Semi-Auto Sheet Press



Thwing-Albert
INSTRUMENT COMPANY

Representante Oficial para América do Sul
REGMED INDÚSTRIA TÉCNICA DE PRECISÃO LTDA

tel: + 55 11 3601-5700

fax: + 55 11 3601-5973

regmed@regmed.com.br

www.regmed.com.br

www.thwingalbert.com