

---

# *IDM*

---

## *Instruments*

### **Thickness Gauge For Foam and Lofty Products**

**Model Number:** F0017

**Details:** This instrument gives an accurate measurement of the thickness to be the basis for accurate values of various properties of cellular materials such as density, tensile strength, tear resistance and compression set. This gauge is suitable for Method A of AS 2282.2.

#### **Applications:**

- Foam
- Lofty products

#### **Features:**

- 80mm circular pressure foot
- 100±10N/m<sub>2</sub> pressure
- Digital scale read-out with an accuracy of ±0.01mm

#### **Options:**

- N/A

#### **Standards:**

- AS 2282.2 Method A

#### **Connections:**

- Not applicable

#### **Dimensions:**

- **H:** 430mm • **W:** 300mm • **D:** 300mm
- **Weight:** 15kg (approx)



**Representante Oficial para América Latina  
REGMED INDÚSTRIA TÉCNICA DE PRECISÃO LTDA.**

**tel. + 55 11 3601-5700 fax. + 55 11 3601-5973 regmed@regmed.com.br www.regmed.com.br**

All sizes and dimensions shown on this data sheet are averages only and images are as accurate as the medium allows.

Products and details shown may change without notice as IDM Instruments reserves the right to make alterations to its product range from time to time.