

## LABORATORY SHEET DRYER TYPE RAPID KÖTHEN *model TD-RK*



The TD-RK is a laboratory sheet dryer designed to be used along with any Rapid Köthen type Laboratory Sheet Former.

This newly designed dryer meets all specification of applicable standards and replaces traditional Rapid-Köthen drying stations saving laboratory space and maintenance costs.

Its elegant design integrates a balanced combination of simplicity, ruggedness and high technology that ensures excellent performance with minimal requirements for installation, operation and maintenance.

### MAJOR TECHNICAL CHARACTERISTICS

|                                 |                                |
|---------------------------------|--------------------------------|
| Drying area                     | 315 cm <sup>2</sup> / Ø 200 mm |
| Drying temperature - adjustable | (90 to 94)°C                   |
| Drying time - adjustable        | (0 to 99) min                  |

*Rugged construction with simple design uses corrosion resistant materials such as stainless steel, bronze and aluminum. Simple and efficient digital controllers for temperature and drying time.*

**Optional accessories:**

· Standard board carriers - Ø 240mm, 250g/m<sup>2</sup>

**Applicable standards:**

· ISO 5269/2



|              |  |
|--------------|--|
| Power supply | 220 VAC, single phase, 1200 W            |
| Air supply   | 2 bar, dry, clean, 0,1 m <sup>3</sup> /h |
| Dimensions   | (65 x 40 x 27) cm (W x L x H)            |
| Weight       | 45 kg                                    |

Note: due to constant development our equipment design and specifications are subject to change without notice.